





General information

- * It is essential that you and any other operator of this product read and understand the content of this manual before installing and using this product !
- * In order to avoid operating errors, these manual must be accessible to the staff at all times.
- * The Hoverboard should only be used for the purpose as described. Please refer particularly to the instructions of the used stretcher.
- * The illustrations do not necessarily correspond to the delivered equipment and are not true to scale.
- * We take no liability for damages caused by operating errors or incorrect assembly or repair.
- * Please pay close attention to the country-related, applicable safety regulations for patient transfer.
- * Subject to technical changes.

Specification

- 1... Pneumatic suspensions with automatic weight adjustment from 40-520 kg/88-1146 lbs and hydraulic shock absorbers.
Total vertical stroke 170 mm/6.7"
- 2... Stretcher mounting platform
- 3... Lateral movement device with pneumatic release (option)
- 4... Switch box, contains:
 - a) 1x Relay for compressor, 12V
 - b) 2x Two-state relays for key circuit, 12V
 - c) 5x Magnetic valves 12V
 - d) Terminal block
- 5... Air pressure tank (approx. 4.5 litres/1.2 gal) with pressure switch
- 6... Reed switch
- 7... Blue momentary main switch
- 8... Red momentary switch for highest (rigid) position (for CPR)
- 9... Momentary switches for pneumatic release of the lateral movement device (option)



Assembly

The Hoverboard may only be installed by qualified personnel (e.g. car technicians or body fitters) and in accordance with the assembly instructions supplied.

The assembly operator is responsible for damages caused by improper assembly !

Technical data

- * Customized Hoverboard for Ferno Viper and Medirol Vivera
 - * Height lowered: 135 mm - 5.3"
Height when active: 220 mm - 9.0"
Height for reanimation (rigid): 310 mm - 12.2"
 - * Standard length: 2003 mm - 78.9"
 - * Total weight 91/119 kg - 200/262 lbs without/with lateral movement device
 - * Maximum loading capacity 520 kg / 1146 lbs (incl. stretcher)
 - * Ignition AND main switch ON: Device ready
Ignition OR main switch OFF: Device lowered (e.g. for loading)
 - * To ensure that the stretcher is always ready for use, the permanent positive is passed through the Hoverboard for continuous charging of the battery.
 - * Electrical connection: **Dimension of all wires must be at least 4,0 mm² (AWG 12) !**
 - Brown = ground (wire no. 31)
 - Red = permanent positive (for charging the Stryker PowerLoad and release of the lateral movement device)
fused in the vehicle with 15A (wire no. 30)
 - Orange = ignition positive, fused in the vehicle with 30A (wire no. 15)
 - Never connect the ignition line together with permanent positive !**
This could lead to consequential damages, which are not covered by warranty !
 - * The valve-control-circuit is internally fused with 5A
 - * Max. power consumption 25A at 12V DC
-

Operation

After successful assembly and electric connection, start ignition of vehicle and turn on the blue main switch.

The Hoverboard now automatically adjusts to the patient's weight and rises to the level for optimal suspension comfort, the switch lights up blue.

The Hoverboard is all-automatic, i.e. when the supply pressure decreases to 8 bar (116 psi), the compressor starts again for approx. 20 seconds, so that there is always sufficient operating pressure.

By switching off the ignition or pressing the main switch again, the Hoverboard lowers automatically for easy loading and unloading. It also lowers by unlocking the stretcher on the Hoverboard.


Pressing the red switch for the CPR-position (for resuscitation) raises the Hoverboard to the highest (fixed) position, the switch lights up red. This function is only active when the main switch is switched on. By pressing it again, the Hoverboard gently lowers to suspended mode.

ATTENTION:

When lifting to the highest position, the compressor runs continuously up to 4 minutes.

Multiple lifting in short intervals will overheat the compressor !


The thermal protection switch will turn off the compressor and it takes approx. 1 hour to cool down !




By pushing one of the pushbuttons, the lateral movement device (option) will unlock pneumatically. While keeping pushed, you can move the Hoverboard crosswise in 8 positions of 32 mm / 1,25" (total = 256 mm / 10").

Releasing the button means locking in the nearest position.

This option also works with ignition off.



Operation



The pneumatic release is supplied by the air tank, which will only be refilled when the main switch is on.
Frequently use while the Hoverboard is switched off leads to temporary malfunction.

In case of malfunction or lack of air pressure you may unlock it by the knobs.

For correct loading and unloading please refer to the instructions of the stretcher.
Please pay attention to the correct locking of the stretcher on the Hoverboard.

Serial number

The serial number is located on the bottom of the plate at the right end and contains the production date.
E.g. serial number **212105** means:
2021 - 21th week, 5th production unit.



Please always quote this number for any complaints or spare parts orders !

Maintenance

One of the many advantages of the Hoverboard-types Airbase, Powerbase, Inbase and Vivibase compared to conventional stretcher tables is that they are not classified as medical devices according to Medical Devices Act (MPG) and are therefore not required to be checked mandatory every year.

Hoverboards are basically maintenance-free, but we recommend an annual inspection with a functional check, regarding e.g. the Powerbase, as part of the maintenance of the stretcher.

This can be managed by one of our certified service partners, whom you find up-to-date on

www.hover.at

You may also contact us directly if you have any questions or problems.

Purification

All Hoverboards are made of high quality stainless steel of the type 1.4301 and are carefully processed in protective atmosphere.

However, „stainless“ does not mean that the material is resistant to all aggressive chemicals, such as e.g. ionized chlorine solutions.

At outside temperatures below -5 ° (23F), calcium chloride is often used as antifreeze solution, which is much more aggressive than the commonly used sodium chloride (“common salt”).

This aggressive solution is brought to the surface of the plate by the wheels of the stretcher and causes surface corrosion as a result.

For a consistently beautiful appearance, it is therefore important to rinse and wipe the surface of the plate daily with clear water during such conditions.

A final impregnation with oil-based stainless steel care products can also help prevent corrosion.

On request, we will be happy to send you a special cleaning and care set.

Disinfection

On the European market, there is such an abundance of disinfectants, solvents and cleaning agents , so that not every single product can be tested.

In addition, the manufacturers are constantly changing and adapting their recipes.

That is why we only use 1.4301 quality for all stainless steel parts, because it means the optimal synthesis of corrosion protection, processability and cost.

The correct (means not too high) concentration of the disinfectant is most important.

Never use disinfectants with chloroacetic acid or other corrosive ingredients !

The bellows is made of polyester fabric with a PVC coating and a PVC support frame. These materials are generally trouble-free, but too high concentrations of disinfectants may cause color damage (bleaching).

The disinfectant also should not remain on the surface anywhere, but should be wiped off after it has acted or also rinsed away with clear water.

Quality

All Hoverboards are tested in accordance with the latest standards by DEKRA in Klettwitz and comply with EN 1789: 2020, EN 1865-5: 2015 and ECE R17 (20 g test).

The Hoverboard Powerbase has also been tested by Stryker Medical EMS (USA) and is homologated for use with the Stryker PowerLoad and PowerPRO XT.

The exclusive use of components of ISO certified pre-suppliers provides industrial manufacturing quality.

Because of our CIP (Continuos improvement programme) and advancement of our products, your Hoverboard may vary from this description.

Spare parts

Air suspension
30620 > with shock absorber
30621 > without shock absorber



30311 Air cushion




30302
Adjustable shock absorber




30057 Magnetic valve block



30058 Magnetic valve 1



30059 Magnetic valve 2



30291 Level control valve



30292 Valve control bow



30272 Air fittings

Spare parts



30211 External compressor
(without Box)

62002 External box
(without compressor)

62003 External box
(with compressor)

30513 JPT-crimp female

30514 JPT-crimp male

30517 JPT-connector 2 pins f

30518 JPT-connector 2 pins m

35060 JPT-connector 6 pins f

35061 JPT-connector 6 pins m



30244 Pressure tube 4 mm, red

30245 Pressure tube 4 mm, blue


30246 Pressure tube 4 mm, yellow

30240 Pressure tube 6 mm, red

30248 Teflon tube,
white

31080 Thermal protection-
tube

30066 Pressure switch,
pre-adjusted 8,5 bar



52014 Switch box, empty


62302 Switchbox for momentary switches,
complete with valves,
terminal block and two-state relays



30072 Reed switch

30073 Block magnet

Spare parts



30083 Bellows, without lateral movement device
30084 Bellows, with lateral movement device



30520 Load relay
30521 Two-state relay
30530 Relay socket



31283 Cable fixer, self-adhesive



35010 Main switch latching (blue)
35013 Main switch momentary (blue)




35011 CPR switch latching (red)
35014 CPR switch momentary (red)




35012 Pushbutton for lateral movement device



30344 Rubber bushing





30361 Plug 60 x 40



35080 Plug combination including contacts



Spare parts



31020 hose connector 4-4-4
 31021 hose connector 4-4
 31022 hose connector 6-6
 31023 hose connector 6-4
 31024 hose connector 8-6

30251 fuse holder
 30252 fuse 30A
 30253 fuse 5A

35022 Connector and cable for front pushbutton
 35023 for rear pushbutton
 35024 for main switch
 35025 for CPR switch



31011 LMD pneumatic cylinder

31030 Linear bearing

31040 Knob

51031 Zahnkamm hinten
 51041 Zahnkamm vorne

51062 Rastbolzen



62260
 Edelstahlpflegeset



51031 Rear toothcomb
 51041 Front toothcomb

51062 Locking bolt

62260
 Care set for stainless steel

Spare parts

Terminal block

- 31270 grey
- 31271 green-yellow
- 31272 blue
- 31274 jumper
- 31275 orange

Terminal block

- 35071 grey
 - 35072 orange
 - 35073 red
 - 35074 blue
 - 35075 black
 - 35076 jumper 2-way
 - 35077 jumper 4-way
-
-



In line with our continuous improvement program, we are always grateful for your suggestions, requests, complaints and proposals for improvement.

Therefore we look forward to hearing from you at

Phone +43 - 660 - 800 9991

Mail info@hover.at

www.hover.at

Have always a safe ride with your new
HOVERBOARD®

Copyright 1/2025

Subject to technical changes

