**AB GERMA®**

## **Vacuum Positioning Pillow**

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# Directions for use

## 2. Precautions

- On a current basis check that the vacuum pump works correctly and that it is possible to get the air out of the Pillow to make it hard.
- Avoid handling the pillow next to sharp objects. If the pillow is pulled be aware of obstacles, sharp edges etc.
- Avoid handling the pillow next to open fire.
- On a current basis check that the pillow has no cracks.
- Own modifications on products not allowed – Your medical specialist to be responsible for any modification.
- When decided to take the product out of service and destroy – make sure to follow your local instructions or contact us for additional information.
- DentaComfort pillow made in compliance with MDR 2017/745 Class 1 concerning Medical Device.

### **Introduction Vacuum Pillow**

- How to use your Germa Vacuum pillow.

The Gernia Vacuum pillow is intended to be used as a tool for, positioning and fixation. Immobilization or stabilisation of patients in many different situations. In order to achieve the desired form and rigidity of the pillow it is necessary to mold it and gradually create vacuum in the pillow by evacuating the air. This can be done by using Germa Vacuum Pump or another kind of vacuum source. The pillow consists of an outer shell of airtight material and a textile inner bag containing granules.

### **Precautions and test**

On a current basis check that the vacuum pump works correctly and that it is possible to evacuate the air from the pillow. Avoid handling of the vacuum pillow on or close to sharp objects or next to open fire. Check for damages, cracks or wear and tear.

### **Specifications**

The vacuum pillow tested and approved according to EN 1021-1. Materials are cover material in PVC-film or PVC coated fabric, inner bag in polyamide and styropor granules. Sizes are: 80 x 100 cm, 80 x 200 cm, 55 x 30 cm, 55 x 65 cm, 55 x 85 cm and 40 x 50 cm and optionally made to customers measurements. Recommended temperature of use -20° C - +40° .

### **Storage**

Recommended storage temperature -20 ° C - +70 °.

### **In use**

Open the valve and put the unit on an even surface. When you have obtained an even distribution of the granules. Shape and fold the unit in its initial position.

Put the vacuum pillow in the desired position on the patient. For pillows with t-valve: Open the valve, evacuate the air and close the valve. For screw valve: Close the valve until it is totally closed – anti clock  
If you are not satisfied with the support and stabilization of the patient. Just open the valve, relocate the pillow in desired position, close the valve, and evacuate again.

### **Cleaning and disinfecting**

After use of the product it can be necessary to rinse it in soap and water. Allow the pillow to dry completely before storing it. If necessary to disinfect. wipe the complete pillow with a solution containing of 70 % alcohol or other most common disinfectants e.g. Lysetol AF.

### **Maintenance and repair**

Test the pillow and pump at regular intervals to confirm the necessary airtightness and functions. Eventual leaks can be found by visual inspection or by reduced filling of air into the pillow followed by

covering it with soapy water all smear the surface and the welds watching out for air bubbles from the pillow. If leaks have been found, clean the area and let it dry. Find a patch in a suitable size and put glue on the patch as well as on the leakage. Wait approx. 5 minutes and place the patch on the leakage. Firmly press the whole area of the patch establishing good contact and press all air out under the patch. The repair will receive full strength after 24 hours.

### **Test**

For testing evacuate the air and create rigidity, close the valve and leave the pillow for about two hours. If it is still as rigid as it was when you evacuated, the pillow may be considered tight. If the pump is not evacuating the pillow, check the pump by closing the tube connector with your fingertip and pump a stroke two-three times. A negative pressure should be created in the tubing and can be felt at your fingertip. Repair kits for pillow is available.

### Accessories and Spare parts

- 21 272 000 701 Germa vacuum hand pump
- 22 230 720 001 Repair kit for pillow

### **Storage**

Store the pillow in room temperature and dry place.

### **Warranty**

Warranty 1 (one) year from delivery date. No liability for outer violence or damages from improper handling.

### **Discard of product**

DentaComfort cushions can be disposed of with normal household waste as specified by local recycling recommendations for plastics according to their active duration. Materials used in the packaging are environmentally friendly and can be recycled.

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**Report no. :**

(To be assigned in by  
Manufacturer)

**Please print or type or fill in by MS Word, all sections a. – i. Fax or e-mail to your supplier.  
Please use repair order when reporting repairs not associated with complaints.**

For claims please specify :  
Purchase Date

**Description of product**

**Description of Incident/Problem:** 1) What, 2) How?

When? In connection with: 1) transport, 2) unpacking, 3) preparing, 4) test, 5) cleaning, 6) normal use, 7) training, 8) disassembly, 9) assembly, 10) storage, 11) other:

**Info on the product, prior to problem,**

Storage: 1) stock, 2) site of use, 3) packed, 4) unpacked, 5) softpack, 6) other:

Used cleaning and disinfecting agents/Sterilizing method:

Pattern of use: 1) stand-by for emergency, 2) routine use, 3) stand by, 4) other:

Environment: 1) Indoor, 2) Outdoor, 3) Heat, 4) Cold, 5) Humidity, 6) Sunlight, 7) Gases, 8) other:

---

Which **actions** have been: a) taken, b) planned by the reporter or user at the site of incident?

- 1) examination of the product on site, 2) replacement of the product, 3) other:

---

Suggested actions by the distributor/complainant

---

---

**Returning of samples/faulty product(s),**

- 1) Already shipped/Will be shipped to \_\_\_\_\_,
- 2) Pending request from manufacturer,
- 3) Product/parts not available:

---

Date:

---

**We acknowledge the receipt of this report by manufacturer and will start processing it immediately. Please use the assigned complaint report number for further reference.**

**Your contact person in this matter is:** \_\_\_\_\_  
Name:

Direct Tel.

Direct Fax.

Init/Date:

---