



# PEDI-MATE PLUS PEDIATRIC RESTRAINT SYSTEM

PediMate<sup>®</sup>Plus Pediatric Restraint System, designed for pediatric patients in sizes from 4.5-45.4 kg. Provides a fast and efficient way to safely secure young patients to almost any cot without having to buy extra equipment or special chairs. PediMate<sup>®</sup>Plus is fully adjustable with a five-point harness system that securely holds patients, which provides comfortable grip during transport. Three straps easily attach PediMate<sup>®</sup>Plus to any crib.

<b>Art. Nr</b>	5003055
<b>Load Capacity</b>	4,5-45 kg
<b>Weight</b>	1,59 kg
<b>Length</b>	90 cm
<b>Width</b>	48 cm



## KEY FEATURES

- Secures pediatric patients in sizes from 4.5-45.4 kg
- Fully adjustable five-point harness system keeps patients safe and provides secure restraint during transport
- Compact rolls for convenient storage and ease of use
- The vinyl construction is non-toxic and easy to clean
- Three restraint straps are easily attached to any crib

