



DECOMPRESSION NEEDLE TPAK 1,6MM

Art. Nr	L845221
Certifications	CE
Color	Red
Weight	14 g
Length	8,2 cm
Width	1,6 mm

Decompression needle for acute treatment of injured people who show symptoms of increased pressure in the alveoli as in hypertensive pneumothorax. The TPAK decompression needle has a length of 82.6 mm, which is longer than the coarsest PVK and therefore penetrates deeper and makes it more likely that you will get it right with thoracic decompression.

TPAK decompression needle is a compact and reliable solution for treating a life-threatening overpressure pneumothorax. TPAK is in a robust packaging with a knurled and easily opened lid. Compact size that is easy to store.



KEY FEATURES

- Easy to use
- 1.6mm (14 gauge) x 82.6mm (3.25")
- Life-saving measure for hypertensive pneumothorax

