



FARETEC INC

RESPONDER PRODUCTS

CT-7 LEG TRACTION SPLINT

ADULT / PEDIATRIC / BELOW KNEE APPLICATION

The CT-7 leg splint has been SIMPLIFIED; applies in 3 easy steps. It has a SHORTER footprint and no longer extends beyond the patient's foot. It has been made more rigid and indestructible with thicker CARBON tubes. The traction mechanism is easy, fast and super precise by means of an acme threaded polycarbonate rod system. We are confident that this leg traction splint is the best choice for any medic.

NSN: 6515-01-706-4198. FDA, MDR/CE. ISO 13485. Made in USA.

Allows stabilization and immobilization of various types of bones fractures in the field. Adjustable in length and circumference. Extends to 1270mm. Max length compressed of 500mm. Stabilizes when compressed. The adjustment/traction point for length and circumference is located at the hip bone when mounted. Weighs 0.6 kg. The splint is radiolucent. Withstands use for over 2 hrs, 3+ times per day for 36 weeks over ten years. Stores at least 10 years in a dry ventilated facility ranging from 0 to 30+C. New CE/MDR according to EU 2017/745 class 1.

IMPROVEMENTS

- No longer extends beyond the patients foot.
- Traction mechanism forces controlled & precise method of traction. Much safer for the patient.
- Simplified layout; no lines, less Velcro, one piece ankle hitch.
- Compact & lighter than the CT-6/any traction splint on the market.
- Increased Diameter of carbon tubing+ 35%. More Rigid fixation.
- Applies traction to patient in under 30 seconds.
- Easier to re-pack, uses new Tactical bag
- One piece ankle hitch with rubber grip/non-slip interior.

SPECIFICATIONS:

SIZE: 12 x 4 x 4in.
WEIGHT: 1 LB
MATERIALS: CARBON TUBING, NYLON & ELASTIC WEBBING
COLORS: BLACK
BAG: TACTICAL BLACK or CAMO
ZIPPER & MOLLE STRAPS



FEATURES & BENEFITS

- ADULTS, PEDIATRICS & BELOW KNEE APPLICATION
 - ACME SCREW TRACTION TECHNOLOGY SUPER PRECISE, STRONG AND ACCURATE.
 - CARBON TUBING. X-RAY SAFE.
 - COMPACT, LIGHT WEIGHT DESIGN.
 - SUPER CLEAN AND EASY TO USE.
- APPLICATION: GOES ON IN 30 SECONDS



1610 West Jackson St.
Unit 6
Painesville OH 44077 USA
+1 440 350 9510
Service@FareTec.com

www.FareTec.com