



LITTER SHIELD



Trusted and proven from years in the field, our Litter Shield attaches to a mount on the litter rail that allows the shield to quickly pivot and nest inside the litter for transporting and storing. A quick-release pin allows the shield to be completely removed when necessary. This lower-profile version makes storage easier, but does not provide as much clearance for the patient's face as the Litter Shield Plus, see below. Each Litter Shield requires a railing clamp.

Replacement parts for Litter Shields are available through the website or by contacting Customer Support.

Option	Litter Shield
Item	740000
PRODUCT WEIGHT	3 LB 2 OZ (1.4 KG)
CLAMP FOR	N/A
Option	19 mm Clamp
Item	740011
PRODUCT WEIGHT	2.7 OZ (77 G)
CLAMP FOR	DISASTER RESPONSE, GAZELLE & JUNKIN LITTERS
Option	22 mm Clamp
Item	740017
PRODUCT WEIGHT	2.4 OZ (68 G)
CLAMP FOR	TRAVERSE ADVANTAGE LITTER
Option	25 mm Clamp
Item	740016
PRODUCT WEIGHT	2.4 OZ (68 G)
CLAMP FOR	CMC STAINLESS & TITANIUM LITTERS
Option	35 mm Clamp
Item	740012
PRODUCT WEIGHT	3 OZ (85 G)
CLAMP FOR	FERNO TYPE, 1-3/8 INCH RAIL
Option	Replacement Pin
Item	740015
PRODUCT WEIGHT	2.1 OZ (60 G)
CLAMP FOR	N/A