



## Features

- Independent leg operation provides safe loading into the ambulance, legs unfold automatically for easy unloading
  - Color-coded controls are easy to use and require minimal training
  - F3<sub>s</sub> Trolley and Stretcher is easy to handle, supports up to 250kg and has a strong structure in aluminum
  - Seven variable heights (controlled by a fully mechanical system) to facilitate loading, unloading and patient transfer
  - Four 200mm diameter swivel wheels with brakes
  - Two front wheels with swivel lock
  - Infinitely – adjustable gas-assist backrest
  - Features raised or contour shock position
  - Swing-down sidearms
  - Four telescopic carrying handles with five positions (longer and shorter included)
  - Four 100mm diameter wheels for the stretcher
  - Foot Support
  - Backrest Airway Management Frame
  - The rugged and lightweight IV Pole offers an integrated height adjustment which allows operators to set their desired height for treatment
  - F3<sub>s</sub> has been designed to comply with and exceed international Standards such as EN 1865-1:2010+A1:2015 and EN 1789:2020
- 60-0350-002 STRETCHER F3S, 5P GRAY
  - 60-0330-004 TROLLEY F3S GRAY

## Accessories & Fasteners



F2 LOCKING DEVICE



MATTRESS  
ORANGE/BLACK



RESTRAINTS KIT



HEADREST



INSTRUMENT TRAY

## Specifications

	Length min.	Length max.	Width	Height min.	Height	Load Height	Weight	Load limit
Stretcher	1945 mm	2350 mm	568 mm		310 mm		18.3 kg	250 kg
Trolley	1890 mm		568 mm	315 mm	850 mm	700 mm	32 kg	250 kg