




## Features

- Independent leg operation provides safe loading into the ambulance, legs unfold automatically for easy unloading
  - Color-coded controls are easy to use and require minimal training
  - F3 Trolley and Stretcher is easy to handle, supports up to 250kg and has a strong structure in aluminum
  - Four 200mm diameter swivel wheels with brakes
  - Two front wheels with swivel lock
  - Features raised or contour shock position
  - Swing-down sidearms
  - Four telescopic carrying handles with five positions (longer and shorter included)
  - Four 100mm diameter wheels for the stretcher
  - Foot Support
  - Secure hardboard (BB) insertion system placed under the cot
  - The rugged and lightweight IV Pole offers an integrated height adjustment which allows operators to set their desired height for treatment
  - F3 has been designed to comply with and exceed international Standards such as EN 1865-1:2010+A1:2015 and EN 1789:2020
- 60-0350-001 STRETCHER F3
  - 60-0330-001 TROLLEY F3 H56
  - 60-0330-002 TROLLEY F3 H50
  - 60-0330-003 TROLLEY F3 H69

## Accessories & Fasteners



F2 LOCKING DEVICE



MATTRESS  
ORANGE/BLACK



RESTRAINTS KIT



HEADREST



INSTRUMENT TRAY

## Specifications

	Length min.	Length max.	Width	Height min.	Height max.	Load Height	Weight	Load limit
Stretcher	1945 mm	2350 mm	568 mm		300 mm		18.3 kg	250 kg
Trolley H56	1890 mm		568 mm	315 mm	850 mm	590 mm	32 kg	250 kg
Trolley H50	1890 mm		568 mm	315 mm	850 mm	540 mm	32 kg	250 kg
Trolley H69	1890 mm		568 mm	315 mm	850 mm	690 mm	32 kg	250 kg