



Features

- Independent leg operation provides safe loading into the ambulance, legs unfold automatically for easy unloading
 - Color-coded controls are easy to use and require minimal training
 - F3S Monobloc is easy to handle, supports up to 250kg and has a strong structure in aluminum
 - Seven variable heights (controlled by a fully mechanical system) to facilitate loading, unloading and patient transfer
 - Four 200mm diameter swivel wheels with brakes
 - Two front wheels with swivel lock
 - The backrest adjusts seven locking position
 - Features raised or contour shock position
 - Swing-down sidearms
 - Four telescopic carrying handles with five positions (longer and shorter included)
 - Foot Support
 - The rugged and lightweight IV Pole offers an integrated height adjustment which allows operators to set their desired height for treatment
 - F3S has been designed to comply with and exceed international Standards such as EN 1865-1:2010+A1:2015 and EN 1789:2020
- 60-0340-004 MONOBLOC F3S, 4P

Accessories & Fasteners



F2 LOCKING DEVICE



MATTRESS
ORANGE/BLACK



RESTRAINTS KIT



HEADREST



INSTRUMENT TRAY

Specifications

Length min.	Length max.	Width	Load Height	Weight	Load limit
1950 mm	2350 mm	568 mm	700 mm	45 kg	250 kg