



## SAFEPAK MULTI-MOUNT READY REFERENCE



SafePak Multi-Mounts attach to iNTRAXX wall tracks to provide an area where small and medium sized SafePak cases may be secured vertically inside the vehicle.

The new Multi-Mounts have been designed with a lower profile and smaller footprint which represents significant weight savings while enabling it to be mounted in more places within the ambulance.

Available in two sizes - 3 Row and 2 Row - and supplied with or without a Safemount connector for the following applications:

### TRACK MOUNTING

- 0480193 3 ROW with Safemount
- 0480201 2 ROW with Safemount

### DIRECT TO WALL MOUNTING

- 0480209 3 ROW plate only
- 0480207 2 ROW plate only

	Multi-Mount 3-Row	Multi-Mount 2-Row
<b>Dimensions</b>	<b>W:</b> 425 mm / 16.7 in <b>D:</b> 40 mm / 1.57 in <i>With Safemount</i> <b>D:</b> 30 mm / 1.2 in <i>Plate only</i> <b>L:</b> 600 mm / 23.6 in	<b>W:</b> 425 mm / 16.7 in <b>D:</b> 40 mm / 1.57 in <i>With Safemount</i> <b>D:</b> 30 mm / 1.2 in <i>Plate only</i> <b>L:</b> 600 mm / 23.6 in
<b>Weight</b>	3.3 kg / 7.3 lb <i>With Safemount</i> 2.5 kg / 5.5 lb <i>Plate only</i>	2.45 kg / 5.4 lb <i>With Safemount</i> 1.6 kg / 3.52 lb <i>Plate only</i>
<b>Load Capacity</b>	12 kg / 26.5 lb	12 kg / 26.5 lb
<b>Crash Tested</b>	26 g	26 g


**SafePaks** are hard-shelled cases lined with pouch mounts which secure SafePouches inside the case. Contents are accessed while mounted to the wall, or removed for use on-scene.

The **Small SafePak** is supplied with a full set of pouches, configured as shown.


The **Medium / Regular SafePak** is supplied with a full set of pouches, configured as shown.



# READY REFERENCE: SAFEPAK MULTI-MOUNT REQUIREMENTS





SafePak Multi-Mount 3 Row



SafePak Multi-Mount 2 Row

Mounting Distance  
350mm between  
hole centres



<b>Part Number</b>	<b>0480024 (Red)</b>
<b>Dimensions</b>	W: 481 mm / 12.8 in D: 240 mm / 9.4 in L: 460 mm / 18.1 in
<b>Weight</b>	7.6 kg / 16.9 lb
<b>Load Capacity</b>	11.3 kg / 25 lb
<b>Crash Tested</b>	10g, 22g, 26g



## Small SafePak

<b>Part Number</b>	<b>0480077 (Blue)</b>
<b>Dimensions</b>	W: 325 mm / 12.8 in D: 240 mm / 9.4 in L: 460 mm / 18.1 in
<b>Weight</b>	5.03 kg / 11.1 lb
<b>Load Capacity</b>	8 kg / 17 lb
<b>Crash Tested</b>	10g, 22g, 26g

Red: 0480023 | White: 0480090

## Medium SafePak

<b>Part Number</b>	<b>0480024 (Red)</b>
<b>Dimensions</b>	W: 481 mm / 12.8 in D: 240 mm / 9.4 in L: 460 mm / 18.1 in
<b>Weight</b>	7.6 kg / 16.9 lb
<b>Load Capacity</b>	11.3 kg / 25 lb
<b>Crash Tested</b>	10g, 22g, 26g

Blue: 0480059 | White: 0480060