



EC Declaration of conformity no.: 140922-001

Product Name:	Restraint System Kangoofix Paediatric Ready2go
Germa Article No:	15271000820
Description:	<p>Kangoofix restraint system intended for use on a child as a harness system for safe transport of the child placed on top of a parent. The parent shall be placed on a stretcher/cot making use of 4-point harness system and upper leg harness where the Kangoofix is fixated. Restraint system consist of five(5) point belt restraints, inner harness of Breath-o-Prene Oeko-Tex 100 with neck/back plate and heat bag with cap.</p> <p>Consist of: 15271000821 Kangoofix 5-point paediatric belt set 15271000822 Kangoofix paediatric Inner harness 15271000823 Kangoofix paediatric Heat bag incl. cap 15271000824 Kangoofix paediatric Neck/Back plate</p> <p>Read User Manual Instructions before use. The system is crash tested and approved according to EN 1789:2007 + A1:2010 and TSFS 2010:2</p>
Accessory:	

Identification: All products with serial numbers issued after Germa LOT number.

We the manufacturer hereby declare that this product is in conformity with the requirements in:

- Annex I and Annex VII and Annex IX in Council directive 93/42/EEC.
Equipment class: I, non-sterile, no measuring function.
- Swedish Medical Devices Act 1993:584

Signed in Kristianstad, Sweden on behalf of AB Germa by:

Position: Quality Manager
Name: Patrik Thornström

Date: 22/09-2014

Sign: 