



EC Declaration of conformity no.: 200629-001

Product Name:	Germa Navy Rescue Stretcher, soft
Germa Article No:	22201811410
Description:	<p>Navy Rescue Stretcher, soft mattress Orange/black, Size 223 x 103x6,5cm, 7 handles on each side + 3 handles at head + 3 handles at foot end</p> <p>Intended use for patient stabilization and handling in difficult environment. It has a head part construction for head stabilization and floating capacity and restraint system. Body part with stabilisation and patient handling and transport of a patient with injuries. Multiple restraint system to secure the patient for evacuation out from a vessel. Navy Rescue Stretcher loading capacity 200 kg. floating support 140 kg</p> <p>NOTE: always a rescuer to close to watch and guard the patient in the equipment. Never left without protection.</p> <p>Including one set of 4 part colour coded restraints with Fix lock buckles and inner protection part with Velcro closing. Rope and strap for navigating and sliding during evacuation</p>
Accessory:	

Identification: All products with serial numbers issued after Germa LOT number.

We the manufacturer hereby declare that this product is in conformity with the requirements in:

- Annex I and Annex VII and Annex IX in Council directive 93/42/EEC.
Equipment class: I, non-sterile, no measuring function.
- Swedish Medical Devices Act 1993:584

Signed in Kristianstad, Sweden on behalf of AB Germa by:

Position: Acting Quality Manager
Name: Björn Holmqvist

Date: 29/6 -2020

Sign: 

Postadress	Adress	Telefon	Telefax	Organisationsnummer. F-skattebevis	Bankgiro	Postgiro
Box 2127 S-291 02 KRISTIANSTAD SWEDEN e-mail: info@germa.se	Industrigatan 54-56 KRISTIANSTAD	+46 (0)44-123030	+46 (0)44-103179	556086-9413	481-4125	21 68 22-7