



VACUUM SPLINT, LONG LEG

Vacuum splints from Ferno[®] are known for their high quality and user-friendly design. The splint versions are among the best-selling worldwide and are manufactured in Sweden, ensuring high standards and affordability.

They are specially designed to stabilize and immobilize fractures in patients during transport. This splint, called full leg, is particularly suitable for immobilizing and stabilizing patients with leg fractures.

| | |
|----------------|----------|
| Art. Nr | GEAS140 |
| Certifications | MDR |
| Weight | 1,15 kg |
| Height | 5,5 cm |
| Length | 112 cm |
| Width | 36/70 cm |



KEY FEATURES

- Made of durable PVC fabric, making them tough and easy to clean
- All seams are fully welded for optimal hygiene
- Splints maintain a vacuum for at least 15 hours, ensuring stabilization over extended periods
- Designed to be radiolucent, allowing scans without removing the splint
- One-year warranty covers manufacturing defects and quality