



THE BAG - 12 PACK - NR3

FOR CHILDREN

The BAG II is a disposable revivator that comes in different sizes.

The transparent mask closes tightly to the patient's face, providing proper ventilation and allowing healthcare professionals to see the patient's lips, secretions and exhaled mist. The joint between the mask and the breathing balloon can be rotated 360 degrees, which gives the healthcare staff the freedom to choose a working position.

The yellow PVC ventilation bladder is transparent and has an ergonomic design. The ergonomic constructions give an excellent feel for how the lungs respond and help to avoid overventilation which can damage the lungs. The BAG II also comes with an adjustable strap that can be attached to the user's hands to reduce physical exertion. Neither the size of the user's hands nor the size of the balloon matters: the breathing balloon attaches to the hand and helps the user to provide controlled ventilation.

Art. Nr	L845221
Certifications	MDD/EC/CE
Color	Yellow/Clear
Liquid Capacity	500 ml
Size	500 ml
Material	PVC, PC, SI, polypropylene

