



FERNO SAVER SAFE EVACUATION CHAIR

Art. Nr	S2001
Load Capacity	180 kg
Weight	11,2 kg
Depth	91,4 cm
Height	134 cm
Width	53 cm
Foldable Height	103 cm

FERNO SAVER SAFE Evacuation Chair provides the ideal solution during an evacuation. The chair is designed for use in emergency situations when it is necessary to safely evacuate persons with reduced mobility down or out of an affected building. The stair chair comes with a cover and mounting brackets.



KEY FEATURES

- Effortless stair descent – reduces manual handling, though operator training is required
- Enhanced control – easy-grip rubber sleeves on handles simplify initial descent positioning
- Patient safety features – adjustable head pad and secure retaining straps
- Stable design – wide rear kickstand adds strength and prevents tipping
- Compact storage – lightweight, with wall bracket, dust cover, and user signage
- Easy maintenance – removable seat cover with underside zips for quick cleaning
- Smooth lateral transfers – seat design supports optimal seated positioning