

Features

Available in Tan

- Fully-Padded Hip Belt
- Includes 1 set of carrying handles as an alternative carrying method.
- Includes 2 adjustment straps for security of product while in transport.

Design

The Traverse Rescue Split Apart Stretcher Backpack protects the basket stretcher from the environmental elements, whether you store it on the top of your vehicle or on a job site. It is designed with ample room to allow for the storage of blankets and other stretcher accessories.

The main zipper is a heavy duty nylon coil zipper with double sliders allowing the top to be completely opened for easy loading and unloading. There is a protective liner and inner padding towards a person's back and on the bottom for personal and equipment protection.

For use with all Traverse Rescue Split Apart stretchers.



Part No.	Description
05-0146	Split Apart Stretcher Backpack, Tan